

# Assessment and differential diagnosis of **PCOS**


## Common symptoms and signs of PCOS

PCOS is defined by a combination of signs and symptoms of androgen excess and ovarian dysfunction in the absence of other specific diagnoses. PCOS signs and symptoms may include the following:

- Hirsutism
- Weight gain
- Acne
- Thinning hair/androgenic alopecia

- Ovulatory dysfunction/menstrual irregularities
- Enlarged polycystic ovaries
- Infertility

## **PCOS Diagnostic Algorithm**



CAH, congenital adrenal hyperplasia; DHEA, dehydroepiandrosterone; FSH, follicle-stimulating hormone; hCG, human chorionic gonadotropin; LH, luteinizing hormone; PCOM, polycystic ovarian morphology; PCOS, polycystic ovary syndrome; TC, test code; TSH, thyroid-stimulating hormone.

PCOM should be assessed by ultrasound. Excluding nonclassic CAH is not indicated for a diagnosis of PCOS phenotype D; testosterone is typically elevated in nonclassic CAH and in range for PCOS phenotype D. If PCOS is suspected, ordering the following with testosterone may facilitate a quicker differential diagnosis: 17-Hydroxyprogesterone (TC 17180), hCG, Total, Quantitative (TC 8396) or hCG, Total with HAMA Treatment (TC 19720)—unnecessary if normal menstrual cycle; Prolactin (TC 746); TSH (TC 899) or TSH with HAMA Treatment (TC 19537); and FSH (TC 470) or FSH and LH (TC 7137). To exclude Cushing syndrome, consider initially ordering cortisol 24-hour urine (TC 14534), saliva (TC 19897), or low-dose Dexamethasone Suppression Test (DST), 1 Specimen (TC 6921).

This figure was developed by Quest Diagnostics based on references below. It is provided for informational purposes only and is not intended as medical advice. Test selection and interpretation, diagnosis, and patient management decisions should be based on the physician's education, clinical expertise, and assessment of the patient.

# Table 1. Laboratory tests for assessment of PCOS criteria and differential diagnosis

Panel components (test code) may be ordered separately.

| Test name (component test codes for panels)   | Test code | Description   |  |  |
|---|-----------|---|--|--|
| Hyperandrogenism  |           |   |  |  |
| Testosterone, Free (Dialysis) and Total, MS<br>Includes total (15983) and free testosterone | 36170     | Diagnose hyperandrogenism   |  |  |
| Testosterone total MS   | 15983     | Diagnose hyperandrogenism   |  |  |
| Differential diagnosis  |           |   |  |  |
| Nonclassic CAH<br>17-Hydroxyprogesterone <sup>a</sup>                                       | 17180     | Diagnose nonclassic CAH   |  |  |
| <b>Hyperprolactinemia</b><br>Prolactin  | 746       | Diagnose hyperprolactinemia                                       |  |  |
| Thyroid disease<br>TSH  | 899       | Diagnose hyperthyroidism or hypothyroidism                        |  |  |
| Amenorrhea caused by hypothalamic amenorrhea or primary ovarian insufficiency               |           |   |  |  |
| FSH (follicle-stimulating hormone)  | 470       | Diagnose primary ovarian insufficiency                            |  |  |
| FSH and LH: includes FSH (470) and LH (615)   | 7137      | Diagnose hypothalamic amenorrhea or primary ovarian insufficiency |  |  |
| Pregnancy<br>hCG, total quantitative  | 8396      | Rule out pregnancy  |  |  |
| Tumors  DHEA Sulfate, Immunoassay <sup>a</sup>  | 402       | Diagnose virilizing tumors  |  |  |

 ${\it CAH, congenital \ adrenal \ hyperplasia; FSH, follicle-stimulating \ hormone; LH, luteinizing \ hormone; PCOS, polycystic ovary \ syndrome; TSH, thyroid-stimulating \ hormone.}$ 

<sup>&</sup>lt;sup>a</sup> This test was developed and its analytical performance characteristics have been determined by Quest Diagnostics. It has not been cleared or approved by the US Food and Drug Administration. This assay has been validated pursuant to the CLIA regulations and is used for clinical purposes.



## Table 2. Polycystic ovary syndrome phenotypes based on Rotterdam Criteria

| Criteria <sup>1-3</sup>                    | Phenotype <sup>4</sup> |          |          |          |  |
|--|------------------------|----------|----------|----------|--|
| Criteria                                   | Α                      | В        | С        | D        |  |
| Hyperandrogenism (clinical or biochemical) | <b>✓</b>               | <b>~</b> | <b>~</b> |          |  |
| Ovulatory dysfunction                      | <b>✓</b>               | <b>~</b> |          | <b>~</b> |  |
| PCOM                                       | <b>✓</b>               |          | <b>~</b> | <b>~</b> |  |

✓ Present

PCOM, polycystic ovarian morphology.



For further information, visit **QuestDiagnostics/PCOS** or scan the QR Code.

#### References

- 1. Joham AE, Norman RJ, Stener-Victorin E, et al. Polycystic ovary syndrome. Lancet Diabetes Endocrinol. 2022;10(9):668-680. doi:10.1016/S2213-8587(22)00163-2
- $2.\,Escobar-Morreale\,HF.\,Polycystic\,ovary\,syndrome:\,definition,\,aetiology,\,diagnosis\,and\,treatment.\,\textit{Nat\,Rev\,Endocrinol}.\,2018;14(5):270-284.\,doi:10.1038/nrendo.2018.244.\,doi:10.1038/nrendo.2018.\,doi:10.1038/nrendo.2018.244.\,doi:10.1038/nrendo.2018.244.\,doi:10.1038/nrendo.2018.\,doi:10.1038/nrendo.$
- 3. Monash University, on behalf of the National Health and Medical Research Council (NHMRC), Centre for Research Excellence in PCOS, and the Australian PCOS Alliance 2018. International evidence-based guideline for the assessment and management of polycystic ovary syndrome 2018. Monash University. Updated October 22, 2018. Accessed August 14, 2022. https://www.monash.edu/\_\_data/assets/pdf\_file/0004/1412644/PCOS\_Evidence-Based-Guidelines\_20181009.pdf
- 4. Azziz R, Carmina E, Dewailly D, et al. The Androgen Excess and PCOS Society criteria for the polycystic ovary syndrome: the complete task force report. Fertil Steril. 2009;91(2):456-488. doi:10.1016/j.fertnstert.2008.06.035

Image content features models and is intended for illustrative purposes only.

## QuestDiagnostics.com

Quest®, Quest Diagnostics®, any associated logos, and all associated Quest Diagnostics registered or unregistered trademarks are the property of Quest Diagnostics. All third-party marks—® and TM— are the property of their respective owners. ©2023-2024 Quest Diagnostics Incorporated. All rights reserved. SB12873 5/2024