

Express Results™

Instant urine drug tests that yield rapid, reliable results

Express Results advantages

- Rapid results, available within a few minutes of the collection
- Accurate and reliable results
- FDA-listed testing options
- Easy to use
- Cup contains built-in specimen validity testing (SVT)
- Variety of drug panel configurations
- Confirmation testing supported at our laboratories

Our proprietary Express Results instant urine testing cup combines the specimen collection and drug test into a single device through the use of an integrated key that begins the screening process once the collection is complete. Results are available for employers within minutes of initiating the test and will remain stable for 60 minutes.

Express Results products are available in a variety of drug panels, and each cup contains built-in specimen validity assurances including testing for oxidizing agents, pH, and creatine, along with an integrated temperature strip. These safeguards help to ensure the integrity of the specimen and the test.

These simple-to-use products are designed to be collected and tested on-site, eliminating the cost and downtime associated with an off-site collection. Express Results kits have an extended shelf life of up to 18 months, and can be ordered directly from Quest Diagnostics once a drug testing account has been established. If the initial drug screen result reads non-negative, the specimen should be sent to one of our laboratories for gas chromatography/mass spectrometry (MS) confirmation testing.

While other instant testing configurations are available upon request, our standard Express Results product line includes the following integrated urine cup testing options and FDA-listed dip cards.



Integrated multi-drug screen cup

Our cup includes all of the features of the multi-drug screen card as well as integrated specimen validity testing. The simplicity of a one-step screening test eliminates the transfer of urine from the collection bottle to the testing cup. Urine is collected in the primary chamber of the device; then the key is activated to separate the urine from the testing chamber to help ensure an uncompromised specimen for confirmation testing. Drug test results are stable for 60 minutes; however, specimen validity testing must be read within three to five minutes of the collection.

- Obtain urine specimen
- Secure lid and push the key to activate testing
- Wait three to five minutes and remove the label on back

Multi-drug screen dip card

Our multi-drug screen dip card tests for 5 drugs and can be used with any collection cup. The dip card yields results in five minutes and delivers accurate results with > 97% correlation at +/-50% of assay cutoffs. It meets urine drug testing requirements for specific industries, offering multiple drug test configurations available. Drug test results are stable for 60 minutes.

- Place card into the urine test specimen
- Wait five minutes to view results
- Results are viewed directly from the test card

Multi-drug screen panels test for the following drugs:

6-panel cup*

Amphetamine

Cocaine

Marijuana

Opiates

PCP

Methamphetamine

10-panel cup*

Amphetamine

Cocaine

Marijuana

Opiates

PCP

Methamphetamine

Barbiturates

Benzodiazepines

Methadone

Propoxyphene

12-panel cup*

Amphetamine

Cocaine

Marijuana

Opiates

PCP

Methamphetamine

Barbiturates

Benzodiazepines

Methadone

Propoxyphene

Oxycodone

Buprenorphine

5-panel dip card

Amphetamine

Cocaine

Marijuana

Opiates

PCP

*All panels include specimen validity testing

For more information about instant drug testing, [visit us at EmployerSolutions.com](http://EmployerSolutions.com).

QuestDiagnostics.com

Quest Diagnostics Incorporated and its subsidiaries (Quest) complies with applicable Federal civil rights laws and does not discriminate on the basis of race, color, national origin, age, disability, or sex. ATTENTION: If you speak English, language assistance services, free of charge, are available to you. Call 1.844.698.1022. ATENCIÓN: si habla Español (Spanish), tiene a su disposición servicios gratuitos de asistencia lingüística. Llame al 1.844.698.1022. 注意: 如果您使用繁體中文 (Chinese), 您可以免費獲得語言援助服務。請致電 1.844.698.1022.

Quest, Quest Diagnostics, any associated logos, and all associated Quest Diagnostics registered or unregistered trademarks are the property of Quest Diagnostics. All third party marks – ® and TM – are the property of their respective owners. ©2018 Quest Diagnostics Incorporated. All rights reserved. SB1562